

SafetyGel® Product Guide

Prevent infection & control fluid waste with
the SafetyGel® super absorbent range



For more information on any SafetyGel® products, contact us on:
T: +44 (0)161 274 1080 | E: sales@safetygel.co.uk | W: www.safetygel.co.uk

Contents

What is SafetyGel®?	3
--------------------------------------	---

Products	4
---------------------------	---

Shaker Pots & Tubs	4
------------------------------	---

Super Absorbent Fast-Acting Sachets	5
---	---

GelPads	6
-------------------	---

Vomit & Urine Bags	7
------------------------------	---

SEE SAFETYGEL IN ACTION

Video Demonstration

Scan to see how rapidly SafetyGel® granules absorb and denature fluid for safe disposal.



OUR SUPER ABSORBENT POLYMER


What Is SafetyGel®?

The SafetyGel® range includes an array of liquid waste control tools. Inside every product are super-absorbent SafetyGel® granules made from our non-toxic formula.

Each gram of SafetyGel® granules is capable of absorbing and denaturing up to 170ml of fluid for simple and safe disposal.

EXAMPLES OF SAFETYGEL CLEANUP SOLUTIONS

Intended Uses




Infection Control
SHARPS CONTAINERS



Hazardous Spill Cleanup
CHEMICALS, BLOOD, URINE



Travel Sickness
LAND, AIR OR SEA TRAVEL



Incontinence Care
BEDPANS, COMMODOES, URINALS

Granules

Shaker Pots ◀ Tubs ▶

SafetyGel® super absorbent & fast-acting granules rapidly convert liquids such as urine, vomit and blood into an easily disposed solid, also absorbing odours.



Shaker Pots & Tubs contain varying volumes of SafetyGel® granules. 1g of SafetyGel® absorbs up to 170ml of fluid.



475g Shaker Pot

SKU: SGS475



1.5kg Tub

SKU: SGTUB15



4.5kg Tub

SKU: SGTUB45

Fast-Acting Sachets


2g Sachets



7g Sachets



SafetyGel® Super-Absorbent Fast-Acting Sachets contain a portion of SafetyGel® granules within our patented water-soluble seal for quick activation in a compact package. Each box contains 100 sachets.



SafetyGel® 2g Sachets, or 'SaferSharps®' are designed for denaturing sharps waste, but are also ideal for small spills from samples.



2g Sachets

SKU: SG2GSAFERSHARPS

Absorption: 340ml



7g Sachets

SKU: SGPS1.2L

Absorption: 1.2L

GelPads

10cm X 10cm



20cm X 20cm



SafetyGel® GelPads rapidly absorb and solidify fluid, trapping bodily waste and odours for secure disposal.

Designed for use with pulp urinals, sick bowls and other clinical waste, GelPads are especially effective for help toileting vulnerable patients or the elderly. Each box contains 100 pads.



10cm X 10cm GelPads

SKU: SGPADVUN1L

Absorption: 1L



20cm X 20cm GelPads

SKU: SGPAD4L100

Absorption: 4L

Vomit and Urine Bags

Vomit & Urine Absorbing Bag

With a capacity of 1000ml, SafetyGel® Vomit and Urine Bags rapidly absorb and trap bodily fluids for secure disposal. Containing the fast-acting SafetyGel® formula, Vomit and Urine Bags also absorb odours.



For after usage, the bags are equipped with an adhesive strip and grip seal to ensure secure containment and disposal.

SKU: SGVOMBAG



Absorb




Contain



Dispose

**Wherever there is a spillage
or hazardous liquid, reach
for SafetyGel® to contain it.**



E: sales@safetygel.co.uk

T: +44 (0)161 274 1080

visit www.safetygel.co.uk